

## Enteric Pathogen Transport Medium

A Cary-Blair modified transport medium suitable for maintenance and enumeration of fecal microorganisms, at room temperature

<b>EPT</b> , 15 ml, 30 ml vials, polypropylene	100/cs	EPT-100
<b>EPT (2 units) &amp; BAG</b> , 15 ml, 30 ml vials, polypropylene	50/cs	EPT-100-2
<b>EPT (3 units) &amp; BAG</b> , 15 ml, 30 ml vials, polypropylene	50/cs	EPT-100-3
<b>EPT</b> , 45 ml, 90 ml container, polypropylene	100/cs	EPT-100C
<b>EPT WITHOUT LABEL</b> , 15 ml, 30 ml vials, polypropylene	100/cs	EPT-101
<b>EPT</b> , 15 ml, 30 ml vials, polypropylene	25/box	EPT-25
<b>EPT &amp; BAG KIT</b> , 15 ml, 30 ml vials, polypropylene	25/box	EPT-25K
<b>EPT/SAF &amp; BAG</b> , 15 ml, 30 ml vials, polypropylene	50/cs	EPT/SAF-50

## SAF Fixative (Sodium Acetate Formalin)

Preserve fecal specimens for all intestinal parasites, exam by wet mount, concentration and permanent stain procedures.

<b>SAF</b> , 15 ml, 30 ml vials, polypropylene	100/cs	SAF-100
<b>SAF</b> , 45 ml, 90 ml container, polypropylene	100/cs	SAF-100C
<b>SAF</b> , 15 ml, 30 ml vials, polypropylene	25/box	SAF-25
<b>SAF &amp; BAG, BULK</b> , 45 ml, 90 ml container, polypropylene	25/box	SAF-25H
<b>SAF &amp; BAG</b> , 45 ml, 90 ml container, polypropylene	25/box	SAF-25T
<b>SAF</b>	4 L bottle	SAF-4L

## Saline 0,85%

Transport of fecal specimens for wet mount exam, especially for motility.

<b>SALINE 0,85%</b> , 15 ml, 30 ml vials, polypropylene	100/cs	SA-100
---	--------	--------

## Formalin 10% Neutral

Preserve fecal specimens for all intestinal parasites, exam by wet mount, and concentration procedures.

<b>FORMALIN</b> , 15 ml, 30 ml vials, polypropylene	100/cs	FN-100
---	--------	--------

## Formalin 10% Buffered

Preserve fecal specimens for all intestinal parasites, exam by wet mount, and concentration procedures.

<b>FORMALIN</b> , 3 ml, 17 X 55 mm	500/cs	FO-03
<b>FORMALIN</b> , 20 ml, 40 ml container, polypropylene	150/cs	FO-20
<b>FORMALIN</b> , 45 ml, 90 ml container, polypropylene	100/cs	FO-45

## PVA Fixative (Poly Vinyl Alcohol)

Preserve fecal specimens for all intestinal parasites and for a number of permanent staining procedures.

<b>PVA</b> , 15 ml, 30 ml vials, polypropylene	PVA-100
<b>PVA WITHOUT LABEL</b> , 15 ml, 30 ml vials, polypropylene	PVA-101

## Vial

Vial, 30 ml, with spoon with or without label, sterile or non sterile, for collection and transport of fresh fecal specimen. 200/bag.

<b>VIAL STERILE</b> , polypropylene	CUC-200S
<b>VIAL</b> , polypropylene	CUC-202

## Kit

Ova and Parasites Concentration System in fecal specimen. 50/cs

**EACH KIT CONTAINS:**     ***Not available in Canada*** PT-50

- Funnel/Filter
- Centrifuge tubes, 15 ml, polypropylene, molded-in grads with plug cap
- Surfactant Triton X-100
- Applicators
- Ethyl Acetate Reagent, 250 ml

**EACH KIT CONTAINS:**     ***Not available in Canada*** PT-50V

- Funnel/Filter
- Centrifuge tubes, 15 ml, polypropylene, printed grads with screw cap
- Surfactant Triton X-100
- Applicator
- Ethyl Acetate Reagent, 250 ml